

1	Status	Submitted
2	Date of Last Update	2019/04/26
3	Person Assigned	David Clunie
4		mailto:dclunie@dclunie.com
5	Submitter Name	David Clunie
6		mailto:dclunie@dclunie.com
7	Submission Date	2019/02/24

8	Correction Number CP-1902	
9	Log Summary: Some IODs that need it are missing the Synchronization Module	
10	Name of Standard	
11	PS3.3	
12	Rationale for Correction:	
13	Some IODs that need it are missing the Synchronization Module.	
14	Correction Wording:	

Amend DICOM PS3.3 as follows (changes to existing text are bold and underlined for additions and ~~struckthrough~~ for removals):

A.53.3 X-Ray 3D Angiographic Image IOD Image Module Table

Table A.53-1. X-Ray 3D Angiographic Image IOD Modules

IE	Module	Reference	Usage
Frame of Reference	Frame of Reference	C.7.4.1	M
	<u>Synchronization</u>	C.7.4.2	<u>U</u>

A.75.3 Parametric Map IOD Module Table

Table A.75-1. Parametric Map IOD Modules

IE	Module	Reference	Usage
Frame of Reference	Frame of Reference	C.7.4.1	M
	<u>Synchronization</u>	C.7.4.2	<u>U</u>