

1	Status	Final Text
2	Date of Last Update	2020/01/17
3	Person Assigned	David Clunie
4		mailto:dclunie@dclunie.com
5	Submitter Name	David Clunie
6		mailto:dclunie@dclunie.com
7	Submission Date	2019/05/09

8	Correction Number CP-1932	
9	Log Summary: Add Series Date and Time to missing Modules	
10	Name of Standard	
11	PS3.3 2019e	
12	Rationale for Correction:	
13	Some Series Modules for IODs that do not the General Series Module do not have the Series Date and Time, e.g., the Encapsulated	
14	Document Series Module.	
15	Correction Wording:	

Amend DICOM PS3.3 as follows (changes to existing text are bold and underlined for additions and ~~struckthrough~~ for removals):

C.7.3.1 General Series Module

Table C.7-5a. General Series Module Attributes

Attribute Name	Tag	Type	Attribute Description
Series Date	(0008,0021)	3	Date the Series started.
Series Time	(0008,0031)	3	Time the Series started.

C.24.1 Encapsulated Document Series Module

Table C.24-1. Encapsulated Document Series Module Attributes

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	The modality appropriate for the encapsulated document. ...
Series Instance UID	(0020,000E)	1	Unique identifier of the Series.
Series Number	(0020,0011)	1	A number that identifies the Series.
<u>Series Date</u>	<u>(0008,0021)</u>	<u>3</u>	<u>Date the Series started.</u>
<u>Series Time</u>	<u>(0008,0031)</u>	<u>3</u>	<u>Time the Series started.</u>
Referenced Performed Procedure Step Sequence	(0008,1111)	3	Uniquely identifies the Performed Procedure Step SOP Instance for which the Series is created. Only a single Item is permitted in this Sequence. Note The Performed Procedure Step referred to by this Attribute is the Step during which this Document is generated.