

THE DICOM 2013 INTERNATIONAL CONFERENCE & SEMINAR

March 14-16

Bangalore, India



Experience of Implementing Radiation Dose Structured Reporting for Fluoroscopy System

Saurav Kumar Baidya

Philips Electronics India Ltd

Software Architect

Bangalore, India

PHILIPS

Experience of Implementing Radiation Dose Structured Reporting for Fluoroscopy System




- **Introduction**
- **Structured Report**
- **Radiation Dose Structured Report**
- **Finding Right Data – A Challenge**
- **Potential Scope of Improvement in Definition**
- **Need of Closer Collaboration**
- **Conclusions**

- **Radiation exposure monitoring is an important aspect of radiation based diagnostic imaging equipments.**
- **This presentation talks about various challenges and observations gathered during the implementation of Radiation Dose Structured Reporting (RDSR) for Fluoroscopy System.**

Structured Report

- Content Items
- Content Items Tree
- Coded content in addition to plain text
- Not presentation guidelines
- Templates



Radiation Dose Structured Report



- **Overcomes weaknesses of MPPS or Image headers as dose monitoring methods**
 - Both MPPS and Image headers lack important dose details
 - MPPS were designed for workflow and are not stored
 - Similar gaps if images deleted or not archived
- **The Radiation Dose Structured Report provide more details in a persistent format**

Radiation Dose Structured Report



TID 10001 (Projection X-ray Radiation Dose)

- * TID 1002 (Observer Context)**

 - # TID 1004 (device Observer Identifying Attributes)**

- * TID 10002 (Accumulated X-ray dose data)**

 - # TID 10004 (Accumulated Projection X-ray Dose)**


- * TID 10003 (Irradiation Event X-ray Data) (one or more)**

 - # TID 1020 (Person Participant) (optional)**

 - # TID 1021 (Device participant)**

Finding Right Data – A Challenge

- **Interpretation of information**
 - Different dose calculations
 - Different positioner primary angles
- **Different sources of data**
 - Patient database
 - Configuration database



Potential Scope of Improvement in Definition

- **Template 10003 row 21 has a 1-n multiplicity**
 - **Rows 22, 23, 24 and 25 are sub elements of 21**
 - **Parent child relationship not clearly defined with respect to the parent with 1-n multiplicity**
 - **Rows 22, 23, 24 and 25 should be combined into a new template**

TID 10003
IRRADIATION EVENT X-RAY DATA
Type: Extensible

21	>	CONTAINS	CONTAINER	EV (113771, DCM, "X-Ray Filters")	1-n	U		
22	>>	CONTAINS	CODE	EV (113772, DCM, "X-Ray Filter Type")	1	U		DCID (10007) X-Ray Filter Types
23	>>	CONTAINS	CODE	EV (113757, DCM, "X-Ray Filter Material")	1	U		DCID (10006) X-Ray Filter Materials
24	>>	CONTAINS	NUM	EV (113758, DCM, "X-Ray Filter Thickness Minimum")	1	U		Units = EV (mm, UCUM, "mm")
25	>>	CONTAINS	NUM	EV (113773, DCM, "X-Ray Filter Thickness Maximum")	1	U		Units = EV (mm, UCUM, "mm")

- **Collaboration with other standards bodies**
- **Spreading radiation dose structured report feature awareness**

- **RDSR implementation is important and experience therefrom is valuable**
- **More collaboration with other bodies will strengthen RDSR acceptance**
- **Need to spread radiation dose structured report feature awareness**
- **Take potential scope of improvement as Correction Proposal**

References

DICOM <http://dicom.nema.org/>

IEC <http://www.iec.ch/>

IHE <http://www.IHE.net/>

Saurav Kumar Baidya

- saurav.baidya@philips.com
- Philips Innovation Campus , Manyata Tech Park, Nagawara Circle, Outer Ring Road , Bangalore – 560045 (Karnataka, India)

Daniel Bauhaus

- daniel.bauhaus@philips.com
- Philips Healthcare, Roentgenstrasse 24, Hamburg, Germany

Mathias Kanitz

- mathias.kanitz@philips.com
- Philips Healthcare, Roentgenstrasse 24, Hamburg, Germany

Thank you for your attention !