

THE DICOM 2013 INTERNATIONAL CONFERENCE & SEMINAR

March 14-16

Bangalore, India



Usage of Presentation State (PR) Objects for Ultrasound 3D images

Ramanjini Dasari

Philips Healthcare

Sr. Technical Specialist
Bangalore, India

Co-Authors: Rupam Sarkar

Software Architect
Bangalore, India

PHILIPS

Agenda


- **DICOM Presentation State**
- **Relationship to Image**
- **3D Presentation State Object**
- **Examples**
- **Conclusions**

- **An independent DICOM Information Object Definition (IOD)**
- **Contains no pixel data**
- **Lightweight object**
- **Saves visual display specification**
- **Multiple presentation states**
- **Manually adjusting visual display settings**

➤ A set of standard DICOM attributes

- (0008,1115) ReferencedSeriesSequence
- (0008,1140) ReferencedImageSequence
- (0008,1150) ReferencedSOPClassUID
- (0008,1155) ReferencedSOPInstanceUID

Relationship to Image (Contd..)




- **Need to transfer studies with 3D images**
- **Implemented DICOM 3D Presentation State object creation**
- **Render the 3D volume same as before capture**
- **Store\Retrieve to\from a PACS system**
- **Use of private tags**

- **3D Presentation State object stores the information given below:**
- **Volume orientation and rotation matrices**
 - **Zoom factor, Gray map, Chroma map (colorization)**
 - **Echo and color vision setting, smoothing, gain**
 - **Brightness, Compression (dynamic range)**
 - **Rotation style (absolute or relative)**
 - **Volume look direction (Top, bottom, left, right, front and back)**
 - **Color baseline, Volume crop info**


3D Presentation State Object (Contd..)

Original image with default Presentation State object




3D Presentation State Object (Contd..)

Same image with different orientation




3D Presentation State Object (Contd..)

Same image with different rotation info




3D Presentation State Object (Contd..)

Same image with different colorization



3D Presentation State Object (Contd..)

Same image with change in Zoom factor



- **An independent DICOM IOD**
- **Helps to display\view DICOM images with different visual display specifications**
 - **3D image volumes in different orientations and rotations**
- **Multiple presentation states**
- **Store\Retrieve to\from a PACS system**
- **DICOM Working Group - Defining an N dimensional Presentation State object**
 - **Enhances interoperability**

References



<http://dicom.nema.org/>

- Ramanjini.Dasari@philips.com
- Rupam.Sarkar@philips.com

Thank you for your attention!